

# Basic procedures for NMR

Dr. E. Manolopoulou, NMR, Lab, Dept of  
Chemistry, Univ. of Crete, March 2019

## Sample Handling

It is a good practice to filter (step 3) NMR solutions directly into the sample tube, or centrifuge them to keep the solution free from dust and other contamination.


**Note:** The sample tube should always be held by the top!

Typical procedures to prepare a sample might be as follows:

1. For a solid sample using a 5 mm tube dissolve up to 20 mg of the sample in about 0.6 ml of the chosen solvent. Typically for a liquid sample, and when observing protons, dissolve 20% sample in 80% deuterated solvent.
2. If the solvent does not contain already, add a small amount (~0.1%) of reference compound Tetramethylsilane (TMS). Make sure the TMS signal is smaller than the most intense sample or solvent signal (otherwise the Signal to noise ratio is wasted because of low receiver gain).
3. If needed, filter the solution into the sample tube through a Pasteur pipette containing a small plug of glass wool.
4. Close the tube with a cap, seal the top with parafilm **to reduce evaporation** and **label the tube** near the top. Be careful to ensure that the cap, parafilm and label are concentric or otherwise they will adversely affect sample spinning

## Sample Preparation for measurement


- Use clean and dry sample tubes (wash tube with EtOH).
- Use medium to high quality sample tubes.
- Always filter the sample solution.
- Always use the same sample volume or solution height.
- Filling volume of a 5 mm tubes is **0.6 ml** or 5 cm.
- Use the sample depth gauge to adjust the sample depth (1.8 cm for older style probes, 2.0 cm for newer style probes).



## Inserting the Sample into the Spinner

1. Sample	4. Center Line
2. Spinner	5. Depth Adjustment Screw
3. Depth Gauge	

- Follow the recommended procedure as shown in the laboratory demonstration.
- Adjust proper position of the sample
- The sample tube should sit tightly inside the spinner.
- Wipe the sample tube clean with EtOH before inserting into magnet or autosampler




**Improper adjustment leads to no detectable signals or breaking the NMR tube inside the detector.**




# Locking the Sample

Lock: for automatic locking

Select solvent



Solvent	Description
Acetic	acetic acid-d4
Acetone	acetone-d6
CS2	benzene-d6
CD2Cl2	dichloromethane-d2
CD3CN	acetonitrile-d3
CD3CN_SPE	LC-SPE Solvent (Acetonitrile)
CD3OD_SPE	LC-SPE Solvent (Methanol-d4)
CDCl3	chloroform-d
CH3OH-D2O	HPLC Solvent (AcetandD2O)
CH3OH-D2O	HPLC Solvent (MethanolD2O)
D2O	deuteriumoxide
D2O_sat	deuteriumoxide with salt
Dioxane	dioxane-d6
DMF	N,N-dimethylformamide-d7
DMSO	dimethylsulfoxide-d6
EtOH	ethanol-d6
H2O+D2O	90%H2O and 10%D2O
H2O+D2O_sat	90%H2O and 10%D2O with salt
H2O+D2O	90%DMSO and 10%DMSO-d6
Juice	fruit juice
MeOD	methanol-d4
Plasma	blood plasma
Py	pyridine-d5
T_H2O+D2O+Me4NO	(CO2)MNO in 90%H2O and 10%D2O, for IMR thermometer
T_H2O+D2O+NaAc	sodium acetate in 90%H2O and 10%D2O, for IMR thermometer
T_H2O+D2O+Pivalate	pivalate-d6 in 90% H2O and 10% D2O, for IMR thermometer
T_MeOD	methanol-d4, for IMR thermometer
TfE	trifluoroethanol-d3
THF	tetrahydrofuran-d3
Tol	toluene-d8
Urine	urine



The lock signal can be controlled **manually** by adjusting Phase, Power, Gain, and Field (with lock off). The Field parameter is Z0.

**Note:** If lock has been achieved, the signal should be of the form of a horizontal line with some associated noise or ripples (see figure below). The height of this line is called the lock level.

# Create a DATA file: Create new experiment

Click or type **edc**

a. **Name:** Define the name to identify your sample. e.g. sample813


b. **EXPNO:** "experiment number" .... Any integer [1 ..... 999].

c. **PROCNO:** 1

d. **Dir:** C:\Bruker\Topspin 1.3\data\username

e. **USER:** Mandatory format "Your-name-month". where month is when the data sets are created.

f. **Solvent:** click and select the solvent in your sample.



**wrpa:** Copy a data set, raw and processed data

**Re:** Read data of specified name or expno

**rpar:** Load experiment parameters

**getprsol:** Get probe and solvent dependent parameters



**NOTE:** Must do getprsol before probe tuning

# Shimming


Shimming is a process in which minor adjustments are made to the magnetic field until the field homogeneity (uniformity) is optimized. Improving the homogeneity will result in better spectral resolution. It will be necessary to reshim each time a sample is changed, especially when the sample solvent changes. The system manager has stored appropriate shim values (in so called shim files) that are frequently updated and will greatly reduce the shimming time required whenever shimming is needed.

Broad lines, asymmetric lines, and a loss of resolution are indications that a magnet needs to be shimmed. The shape of an NMR line is a good indication of which shim is misadjusted. Consider a single narrow NMR line. If we zoom in on this line we might see the following shape. The following series of spectra depict the appearance of this spectral line in the presence of various inhomogeneities.

rsh: read a recent shimfile

wsh: write shims


## Shimming -errors



*You can check the line shape with a quick single scan to obtain a <sup>1</sup>H spectrum after any shimming*

# Manual Shim

1. Take note of the lock level in the lock display as the initial reference level.
2. Start with the functions **Z** and **Z<sup>2</sup>**, **Z<sup>3</sup>**
3. Readjust **Z** and **Z<sup>2</sup>**
4. For the first sample that you do, adjust also functions **X** and **Y** then **XZ** and **YZ** (if there is a large change return to X and Y), then **XZ<sup>2</sup>** and **YZ<sup>2</sup>** (if there is a large change return to **X, Y, XZ** and **YZ**) and then **XY** and **X<sup>2</sup>-Y<sup>2</sup>**
5. Then you may adjust **Z<sup>3</sup>** and return to adjust **Z** and **Z<sup>2</sup>**
6. Finally, collect a quick one scan proton spectrum and review for sharp and symmetrical line shape



For best resolution, adjust the X, Y and XY shims to obtain the maximum lock level.

For most simple samples, steps 1-3 will be sufficient to obtain a good spectrum



## Tuning and Matching the Probe

**Tuning** involves adjusting the probe circuitry so that the frequency at which it is most sensitive is the relevant transmission frequency (SFO1, SFO2 etc.) Each coil in the probe will be tuned (and matched) separately. If the probe has been changed or the transmission frequency altered significantly, it may be necessary to retune the probe. For routine work in organic solvents with selective probes, the value of the transmitted frequencies are unlikely to vary greatly

### Matching


involves ensuring that the maximum amount of the power arriving at the probe base is transmitted up to the coil which lies towards the top of the probe. This ensures that the minimum amount of the power arriving at the probe base is reflected back towards the amplifiers (and consequently wasted).

**Note:** Bruker offers two different types of tuning and matching adjustments


**atma:** for probe with automatic tuning THIS is used on the 300 MHz

**atmm:** Tune/match ATM probe manually

**atma exact:** will determine the optimum tuning and matching more precisely than atma without an argument and will therefore be slower




# Examples of Wobble Curves with Different Tuning and Matching



1.	Bad matching and tuning.	3.	Good matching, bad tuning.
2.	Bad matching, good tuning.	4.	Good matching and tuning.

## Receiver Gain

The receiver gain is a very important parameter that is used to match the amplitude of the FID to the dynamic range of the digitizer. The gain is set by clicking **rg** and providing a suitable value



The NMR-signal received from the resonant circuit in the probehead needs to be amplified to a certain level before it can be handled by the computer.

The detected NMR-signals vary over a great range due to differences in the inherent sensitivity of the nucleus and the concentration of the sample.

Rg values: 0-200

**Rga**: automatic receiver gain optimization

**In all cases, the automatic receiver gain command should be used on the 300 MHz**


# Acquisition parameters

**d1** – relaxation delay for proper integration.

**td** –the number of points to be sampled and digitized to form the FID.

**Final number of accumulated scans = NS +DS**

Typically TD is set to 16, 32 or 64K for standard 1D FID's. Although increasing the size of TD will improve the FID resolution, it does so at the expense of a longer acquisition time.



**ns**: number of scans. Increasing the number of scans increases signal-to-noise ratio, but also increases overall experiment time.

**ds**-number of dummy scans. Several sets of pulses which are identical to those used for acquisition are sometimes transmitted to the sample before any signals emitted by the sample are allowed to enter the receiver. **This is to allow the sample to reach a steady state or equilibrium.** The number of such dummy scans depends on the sample relaxation time and susceptibility to heating. Typical values for standard experiments are 4 or 8. The pulse program can be consulted in other cases.


**o1**: Transmitter Frequency **Offset** for Channel F1 in **Hz**. This is the frequency used to excite the observe nucleus (logical channel F1) and will be at the center of the spectrum. It may be thought of as the central frequency in the window through which the spectrum is observed

**o1p**: Transmitter frequency offset for channel F1 in **ppm**

# Check your acquisition parameters

**used** - check and optimize experimental setup (Display important parameters)

Modify acquisition parameters if desired.



Pulse Program. zg30 means 30° hard pulse, and acquire

Number of points collected. Usually 32k or 64k

Sweep Width, in Hz and in ppm

Relaxation Delay, usually 2 seconds for

Number of scans

Carrier Frequency in ppm. Center of your spectrum


90° pulse length and power. The **getprosol** command will populate these values for you.

You can also use the command **pulsecal**, which will find the accurate 90° pulse

**eda** - To access all of the acquisition parameters

**Note:** The disadvantage of eda is that it shows many parameters, most of which do not have to be changed for a standard experiment. Generally, it is more convenient to use used to set the acquisition parameters. This command only shows the parameters which are actually used for the current experiment

Click the AcqPars tab to display the acquisition parameters




 show the pulse program parameters only

 button to read the prosol parameters.

To start the acquisition:

– Click  in the upper toolbar or enter **zg on the command line.**

The data window toolbar will automatically switch to the Acqu tab and the FID display window will appear:




# Acquire data

**zg-**record a spectrum

**multizg.-**(for # experiments) It will ask you for the number of experiments to be performed. If this number is higher than the number of experiments you have set up, it will repeat the last one for the remaining experiments.

**expt:** display the experiment time



The buttons in the toolbar have the following functions:

-  Show FID in shuffled mode.
-  Show FID in unshuffled mode, horizontally arranged.
-  Show FID in unshuffled mode, vertically arranged.
-  Show FID in unshuffled mode, interleaved.
-  Switch between FID and spectrum.
-  Stop the acquisition [ **stop**].
-  Halt the acquisition [ **halt**].




**tr- optional step.** to preview an on-going experiment (such as a C-13), before its completion

**halt:** To halt the measurement *before it completes the NS scans*

**DO NOT TYPE STOP!** **Stop** command stops the acquisition immediately without writing the data to disk


**Kill:** kill a specific process

# Processing 1D data



Action	Command	Comment
exponential multiply	<b>em</b>	lb is linebroadening in Hz
fourier transform 1D	<b>ft</b> <b>fp</b> <b>ef</b> <b>efp</b>	fourier transform only ft+pk em+ft em+ft+pk
Automated phase correction	<b>apk</b>	Not always sufficient
Baseline correction	<b>abs</b>	Automatic baseline correction and integration on a 1D-spectrum
	<b>absn</b>	Add the letter 'n' after the command to avoid the automatic setting of integral limits.
Automatically calibrate the spectrum	<b>sref</b>	calibrates the TMS peak to 0 ppm


# 1D Interactive Window Multiplication Procedure



You can perform interactive window manipulation as follows:

1. Select the window function (parameter WDW).
2. Set the corresponding parameter(s), e.g.
  - LB for exponential.
  - LB and GB for Gaussian.
  - SSB for sine bell and squared sine.

The displayed spectrum and/or FID will be automatically adjusted as you change the window function and parameters.




3. Click the Save button to store the window settings and return. Now you can perform further processing steps like Fourier transform, phase correction etc.

Exponential multiplication (**em**): can improve S/N. The default lb is 0.3 Hz; employing a higher line broadening value can yield better S/N, but at the expense of resolution.

Gaussian multiplication (**gm**): used for resolution enhancement using lb as a negative value e.g. -1, -2 Hz and gb (Gaussian broadening; 0 to 1); typically 0.2-0.3.

This reduces linewidths, making multiplets clearer, but at the expense of S/N and may also lead to peak-shape artifacts

examples



# How to use the Toolbar (1D)

The lower toolbar contains buttons for display functions.

## Buttons for vertical scaling (intensity)




**\*2** Increase the intensity by a factor of 2 [**\*2**].

**/2** Decrease the intensity by a factor of 2 [**/2**].

**↕** Change the intensity smoothly.

## Buttons for interactive functions



The functions of the individual buttons are:

**↶** Enter phase correction mode.


**↷** Enter calibration mode.

**⤿** Enter baseline correction mode.

**⊥** Enter peak picking mode.

**∫** Enter integration mode.

**⊞** Enter multiple display mode.



## Buttons for horizontal shifting



**←** Shift to the left, half of the displayed region [**.sl**].

**↔** Smoothly shift to the left or to the right.

**→** Shift to the right, half of the displayed region [**.sr**].

**⏪** Shift to the extreme left edge of the spectrum [**.sl0**].

**⏩** Shift to the extreme right edge of the spectrum [**.sr0**].

## Buttons for vertical shifting



**⏴** Shift the spectrum baseline to the middle of the data field [**.su**].

**↕** Smoothly shift the spectrum baseline up or down.


**⏵** Shift the spectrum baseline to the bottom of the data field [**.sd**].

# Phase correction

**pk**: apply phase correction (apply **phc0** and **phc1** to spectrum)

**apk**: automated phase correction

## Manual phase correction



Right click on the most upfield or most downfield peak, and select Set Pivot Point.


To phase, click and hold the 0 or the 1 order phase, and move up or down. Use the 0 order for the peaks near the pivot point, and use the 1 order for the remainder of the peaks.

When you are done, hit the save and close button.

# Peak picking

pp- Perform peak picking


Or select the Peak Picking button under the Process tab . The peak picking window will open.



cancels the selected peaks




allows individual peak picking



# How to Perform Interactive S/N Calculation

The current signal region (parameters SIGF1-SIGF2) and noise region (parameters NOISF1-NOISF2) are displayed.



Data window in S/N measurement mode.

2. Move the mouse into the data window.  
3. Left-click-hold and drag the mouse from one edge of the *signal region to the other edge*.

A horizontal double-headed arrow will indicate the signal region.

4. Left-click-hold and drag the mouse from one edge of the *noise region to the other edge*.

A horizontal double-headed arrow will indicate the noise region.


5. Right-click any position in the data window. The popup menu as shown will appear

Quit
Clear NOISEREG
Clear SIGREG
Edit Regions ...
Change Region Widths ...
Start S/N Calculation
Enter Zoom

Choose **Start S/N calculation**.



The other entries allow you to redefine or clear the regions. After the S/N calculation has finished, the result will appear on the screen.

# Integration




You can integrate 1D spectra by selecting manual integrate from Analyse tab/integrate

From here you can see all your integrals




# Incorrect Integration


Correct phase and baseline



Incorrect baseline (sloping)






Incorrect phasing



# Processing 2D data (auto-commands)

For COSY



**xfb**: process data including FT in both directions F1 and F2

or **xf1**: process data including FT in F1 and

**xf2**: process data including FT in F2



**sym**: symmetrise spectrum about the diagonal

**abs1**: automatic baseline correction in F1

**abs2**: automatic baseline correction in F2

# Processing 2D data (auto-commands)

For HSQC & HMBC




**xfb**: process data including FT in both directions F1 and F2  
or **xf1**: process data including FT in F1 and  
**xf2**: process data including FT in F2

**abs1**: automatic baseline correction in F1  
**abs2**: automatic baseline correction in F2  
**apk2d**: automatic 2d phasing



# Processing 2D data

## Advanced 2D NMR data processing with TOPSPIN



Manual phase correction

# Manual phase correction




Select manual phase correction from process tab/adjust phase


Place mouse to one of the major signals, right click and "add"

Repeat the same procedure to select another signal far away from the previous one, right click and "add"

Click  in the new window



# 2d phase correction on rows



Use 0 and 1 for phase adjustment



Click save when finished with rows phase adjustment

Repeat same phasing procedure for columns by clicking



After phasing the column click save and return

Finalize the phase correction



# Bibliography

NMR Spectroscopy (Some technical and practical aspects)

University of Oxford  
Doctoral Training Centre



Measuring Proton NMR Spectra by Roy Hoffman and Yair Ozery 20th September 2010

User Guide for Bruker NMR Spectrometers

TopSpin User Manual Version 003 by Bruker

Basic NMR Hands-on Operation Manual by U. of Maryland, NMR Center, Fu Chen

Processing With Topspin by Duncan Howe October 31, 2016 ,Department Of Chemistry  
University Of Cambridge

<http://triton.iqfr.csic.es/guide/nmr/manual>

[www.uleth.ca/artsci/sites/artsci/files/docs/2D-TrainingManual.pdf](http://www.uleth.ca/artsci/sites/artsci/files/docs/2D-TrainingManual.pdf)

UC Davis NMR Facility Bruker short manual for Topspin 3.X

A User Guide to Modern NMR Experiments by Barbara Odell and Tim Claridge, University of Oxford.