


# Basic 1D experiments

Dr. E. Manolopoulou, NMR, Lab, Dept of  
Chemistry, Univ. of Crete, March 2019

# 1D <sup>1</sup>H NMR

## PRELIMINARY SET-UP

1. Insert the sample to the spectrometer
2. Choose the right deuterated solvent with lock command
3. Create a new dataset (**edc**) and read the standard BRUKER parameter set (**rpar**) to record a conventional <sup>1</sup>H spectrum with **rpar PROTON all** (the pulse program **zg30** can be visualized in the **PulsProg** section or with the **edcpul** command).
4. **getprsol**: Get probe and solvent dependent parameters (corresponding pulses and power levels)
5. Tune and match the probehead (**atma** or **atma exact**)
6. Optimize the shim procedure (read an optimized shim file with the **rsh** command and perform shimming)



By default, the following parameters are set to:


Relaxation delay (1-5\*T1) (**d1**): 1s  
Number of scans (**ns**): 16  
Dummy scans (**ds**): 2  
Spectral Width (**sw** in ppm): 20.8  
Center of spectrum (**o1p** in ppm): 6.175  
Time Domain (**td** in ppm): 32k

## ACQUISITION

7. Set the appropriate **ns** and **ds** for the experiment
8. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

## PROCESSING


9. The recorded data is Fourier transformed with **ft** (or **ef**) and phase and baseline corrections are performed using **apk** and **absn**, respectively.



# 1D <sup>1</sup>H NMR with water presaturation



## Step 1: PRELIMINARY SET-UP

1. Insert the sample to the spectrometer
2. Choose the right deuterated solvent with lock command
3. Create a new dataset (**edc**) and read the standard BRUKER parameter set (**rpar**) to record a conventional <sup>1</sup>H spectrum with **rpar PROTON all** (the pulse program **zg30** can be visualized in the **PulsProg** section or with the **edcpul** command).
4. Tune and match the probehead (**atma** or **atma exact**)
5. Optimize the shim procedure (read an optimized shim file with the **rsh** command and perform shimming)
6. Record a typical <sup>1</sup>H spectrum. Note the frequency (**o1**) of the solvent resonance



## Step 2: Put the solvent peak on-resonance by:

1. expand about the solvent peak enough that you can easily see the center
2. click  and then left-click with the cursor in the middle of the solvent peak
3. choose o1
4. Write down the value for o1 in Hz.




## SPECIFIC PARAMETERS (acqupars)

- The power level and the duration of the presaturation are defined by
- **p19** (start with **55 dB**) and **d1** (**2s**), respectively.
- minimum number of scans **ns=8 ds=4**
- **p19** (**38-40 dB**)
- **Td** (**32k**)

**Step 3:** 1. Create a new dataset with **edc** and change the pulse program (**pulprog zgpr**).

2. On the command line, type "o1" and key in the solvent frequency that you get from the previous experiment

3. Type **rga**
4. Set **ns, ds**
5. Type **zg** to start acquisition
6. Process spectrum (**ef,apk,absn**)



# The Standard $^{13}\text{C}$ NMR Experiment

## Zgdc30

1D-sequence **with decoupling**, using a  $30^\circ$  flip angle. Result is a standard  $^{13}\text{C}$  NMR spectrum with proton broad-band decoupling

### PRELIMINARY SET-UP


1. Insert the sample.
2. Choose the solvent deuterium signal with the **lock** command.
3. Check shimming from  $^1\text{H}$  spectra
4. Create a new dataset with **edc** and read the standard parameter set to record a conventional  $^{13}\text{C}$  spectra with **rpar C13CPD32 all** (the pulse program **zgdc30** can be displayed with the command **edcpul**)
5. Update the corresponding pulses and power levels in the acquisition parameters according to the selected solvent/probehead parameters by executing the **getprosol** command
6. Tune and match the probehead (**atma**)

### ACQUISITION

7. Set the appropriate **ns** and **ds** for the experiment
8. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

### PROCESSING

9. The recorded data is Fourier transformed with **ef** and (**lb=1**) and phase and baseline corrections are performed using **apk** and **absn**, respectively.



### SPECIFIC PARAMETERS (acqupars)


By default, the following parameters are set to:

**ns=32, ds=2**

**d1=2**

**sw=331.2**

**td=64k**




# Quantitative measurements with inverse gated methodology

## 1D <sup>1</sup>H-decoupled <sup>13</sup>C spectrum without NOE

### PRELIMINARY SET-UP

1. Insert the sample.
2. Choose the solvent deuterium signal with the **lock** command.
3. Check shimming from <sup>1</sup>H spectra.
4. Create a new dataset with **edc** and read the standard parameter set to record a <sup>1</sup>H-decoupled <sup>13</sup>C spectrum without NOE with **rpar C13IG all** (the pulse program **zlg30** can be displayed with the command **edcpul**)
5. Update the corresponding pulses and power levels in the acquisition parameters according to the selected solvent/probehead parameters by executing the **getprosol** command
6. Tune and match the probehead (**atma**)



### SPECIFIC PARAMETERS (acqparams)

By default, the following parameters are set to:


**ns**=16 **ds**=4,  
**d1**=60s  
**sw**=250  
**td**=32k

### ACQUISITION

7. Set the appropriate **ns** and **ds** for the experiment
8. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

### PROCESSING

9. The recorded data is Fourier transformed with **ef** and (**lb**=1) and phase and baseline corrections are performed using **apk** and **absn**, respectively.




# DEPT -135 <sup>13</sup>C NMR experiment

## Full decoupled <sup>13</sup>C spectrum

### PRELIMINARY SET-UP

1. Insert the sample.
2. Choose the solvent deuterium signal with the **lock** command.
3. Record a conventional <sup>1</sup>H-decoupled <sup>13</sup>C spectrum.
4. Create a new dataset with **edc** and read the standard parameter set to record a DEPT spectrum with **rpar C13DEPT135 all** (the pulse program **dept135** can be displayed with the command **edcpul**).
5. Update the corresponding pulses and power levels in the acquisition parameters according to the selected solvent/probehead parameters by executing the **getprosol** command
6. Tune and match the probehead (**atma**)



### SPECIFIC PARAMETERS (acqupars)

By default, the following parameters are set to:

**ns**=128, **ds**=2

**d1**=2s, **d2**=3.57ms

**sw**=331.2


**td**=64k

### ACQUISITION

10. Set the appropriate **ns** and **ds** for the experiment
11. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

### PROCESSING


12. The recorded data is Fourier transformed with **ef** and (**lb**=1) and baseline correction **absn**.



# 1D $^{31}\text{P}$ -decoupled $^1\text{H}$ spectrum

## PRELIMINARY SET-UP

1. Insert the sample in the spectrometer
2. Create a new dataset with **edc** and read the standard parameter set to record a DEPT spectrum with **rpar PROP31DEC all** (the pulse program **zgig30** can be displayed with the command **edcpul**).
3. Update the corresponding pulses and power levels in the acquisition parameters according to the selected solvent/probehead parameters by executing the **getprosol** command



## SPECIFIC PARAMETERS (acqupars)

By default, the following parameters are set to:

**ns**=16 **ds**=2,

**d1**=1s

**Sw (ppm)** =150

**td**=32k


**o1p (ppm)**= 6.175, **o2p**= center of the  $^{31}\text{P}$  spectrum

## ACQUISITION

5. Set the appropriate **ns** and **ds** for the experiment
6. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

## PROCESSING

7. The recorded data is Fourier transformed with **ef**, phase and baseline corrections are performed using **apk** and **absn**, respectively.






The  $^{31}\text{P}$  chemical shift range is rather large and covers approximately from -180 to 250. The default sweep width of the Bruker standard  $^{31}\text{P}$  parameter sets may not cover the whole chemical shift range and adjustment may be needed.

# Selective 1D Experiment: NOE

## PRELIMINARY SET-UP

1. Run a **1D Proton spectrum**
2. Process with **efp**, **apk**.
3. Identify your target for your 1D NOESY
4. Define your regions:
  - a. Under Process tab, hit the Integration button.
  - b. Delete all integrals if any are present
  - c. With the integration cursor enabled, click and drag over your peak to define the region.
  - d. Hit save as, save to region, then save and close



## SPECIFIC PARAMETERS (acqupars)

By default, the following parameters are set to:

**ns**= 32, 64, 128 **ds**=?

**d1** =3s, **d8** (NOESY mixing time)= 0.1-0.8s (for large molecules to small ones)

**sw**=?

**td**=?

5. Create a new dataset with **edc** and read the standard parameter set to record a 1D NOESY spectrum with **rpar SELNOGP**
6. Update the corresponding pulses and power levels in the acquisition parameters according to the selected solvent/probehead parameters by executing the **getprosol** command



## ACQUISITION

7. Set the appropriate **ns** and **ds** for the experiment
8. In the sel1d message window, click **OK to start the acquisition**

## PROCESSING

9. The recorded data is Fourier transformed with **ef** and phase and baseline corrections are performed using **apk** and **absn**, respectively.

## 1D NOESY



## Bibliography

1D and 2D Experiments Step-by-Step Tutorial. Basic Experiments User Guide (Version 004) Bruker

<https://pharm.ucsf.edu/sites/pharm.ucsf.edu/files/AVANCE%20Beginner%27s%20Guide.pdf>

TopSpin Guide Book. Basic NMR Experiments User Manual (Version 002) Bruker

[https://www.nmr.ucdavis.edu/sites/g/files/dgvnsk4156/files/inline-files/iconnmruserguide-ucdavis\\_chem.pdf](https://www.nmr.ucdavis.edu/sites/g/files/dgvnsk4156/files/inline-files/iconnmruserguide-ucdavis_chem.pdf)

1D and 2D Experiments Step-by-Step Tutorial Advanced Experiments User Guide (Version 002) Bruker

<http://triton.iqfr.csic.es/guide/tutorials/multnuc/h1dec31.html>