



Basic 2D NMR experiments

Dr. E. Manolopoulou, NMR, Lab, Dept of
Chemistry, Univ. of Crete, March 2019

^1H - ^1H COSY

PRELIMINARY SET-UP

1. Insert the sample to the spectrometer
2. Choose the right deuterated solvent with lock command
3. Create a new dataset (**edc**) and read the standard BRUKER parameter set (**rpar**) to record a conventional ^1H spectrum with **rpar PROTON all** (the pulse program **zg30** can be visualized in the **PulsProg** section or with the **edcpul** command).
4. **getprsol**: Get probe and solvent dependent parameters (corresponding pulses and power levels)
5. Tune and match the probehead (**atma** or **atma exact**)
6. Optimize the shim procedure (read an optimized shim file with the **rsh** command and perform shimming)
7. Record a typical ^1H spectrum. Note the **SW** values to optimize spectral widths in the corresponding 2D experiment.
8. Create a new dataset with **edc** (or read and copy a pre-existing experiment) and read the standard BRUKER parameter set (**rpar**) to record a 2D ^1H - ^1H COSY-45 spectrum with **rpar COSYGPSW all** (the pulse program **cosygpqf** can be visualized in the **PulsProg** section or with the **edcpul** command)



SPECIFIC PARAMETERS (acqparms)


By default, the following parameters are set to:
2 td in F2 (1K-2K), **1 td** in F1 (128w- 256w)
ns=8, ds=4
1 sw= 12, 2 sw =12
d1=2

ACQUISITION

11. Set the appropriate **ns** and **ds** for the experiment
12. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

PROCESSING


13. Process the recorded data with **xfb**
14. The resulting 2D spectrum can be symmetrized by using the **sym** command.




9. **getprsol**: Get probe and solvent dependent parameters (corresponding pulses and power levels for the current experiment)
10. COSY experiments must be run without sample spinning

^1H - ^1H COSY with presaturation

PRELIMINARY SET-UP

1. Perform steps **1-7** as previous
2. Find the ω_1 of the solvent resonance
- 3. Put the solvent peak on-resonance by:**
 - a. expand about the solvent peak enough that you can easily see the center
 - b. click  and then left-click with the cursor in the middle of the solvent peak
 - c. choose ω_1
 - d. Write down the value for ω_1 in Hz.



4. Create a new dataset (**edc**) (or read and copy a pre-existing experiment) and read the standard BRUKER parameter set (**rpar**) to record a 2D ^1H - ^1H COSY spectrum using presaturation change the pulse program to **cosygpprqf** that can be visualized in the **PulsProg** section or with the **edcpul** command).

5. **getprsol**: Get probe and solvent dependent parameters (corresponding pulses and power levels for the current experiment)

On the command line, type " **ω_1** " in the solvent frequency that you get from the previous experiment

SPECIFIC PARAMETERS (acqupars)

By default, the following parameters are set to:

2 td in F2 (1K-2K), **1 td** in F1 (128w- 256w)

ns=1, ds=8

1 sw= 12, 2 sw =12

d1=1.5-2s


pl9=55-60dB

ACQUISITION

7. Set the appropriate **ns** and **ds** for the experiment
8. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

PROCESSING

9. Process the recorded data with **xfb**
10. The resulting 2D spectrum can be symmetrized by using the **sym** command.




^1H - ^{13}C HSQC phase-sensitive

HSQCETGPSI (hsqcetgpsi) – simple gradient HSQC, non-Edited, sensitivity improved (si).

PRELIMINARY SET-UP


1. Run a **conventional ^1H spectrum**. Note the **o1p** and **SW** values to optimize spectral widths in the corresponding 2D experiment
2. If required, record a ^1H -decoupled ^{13}C spectrum
3. Create a new dataset (**edc**) (or read and copy a pre-existing experiment and skip step 2) and read the standard BRUKER parameter set to record a phase-sensitive ge-2D ^1H - ^{13}C HSQC spectrum. Change the pulse program to **hsqcetgpsi** can be visualized in the **PulsProg** section or with the **edcpul** command).



4. **getprsol**: Get probe and solvent dependent parameters (corresponding pulses and power levels for the current experiment)
5. HSQC experiments must be run without sample spinning
6. Tune and match the probehead (**atma**)

ACQUISITION

7. Set the appropriate **ns** and **ds** for the experiment
8. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).




SPECIFIC PARAMETERS (acqupars)

By default, the following parameters are set to:

- 2 td** in F2 (1K-2K), **1 td** in F1 (64w-256w)
- ns=96, ds=16**
- 1 sw**= 250 ppm (F1), **2 sw** =12 ppm (F2)
- d1=1.5-2s**

PROCESSING

9. Process the recorded data with **xfb**
10. The resulting 2D spectrum can phase and baseline corrected by **apk2d** and **abs1, abs2**




1H-13C long range coupling HMBC

HMBCGP (hmbcgpndqf) – Gradients for coherence selection (gp), low pass filter (lp), no decoupling during acquisition (nd), and magnitude mode (qf). Simple and no 180° pulses.

PRELIMINARY SET-UP

1. Run a **conventional 1H spectrum**. Note the **o1p** and **SW** values to optimize spectral widths in the corresponding 2D experiment
2. If required, record a 1H-decoupled 13C spectrum (**o2p, SW**)
3. Create a new dataset (**edc**) (or read and copy a pre-existing experiment and skip step 2) and change the pulse program to **hmbcgpndqf** in the **PulsProg** section or with the **edcoul** command.




SPECIFIC PARAMETERS (acqupars)

By default, the following parameters are set to:

- 2 td** in F2 (1K-2K), **1 td** in F1 (64w-256w)
- ns=2, ds=16**
- 1 sw**= 250 ppm (F1), **2 sw** =12 ppm (F2)
- d1**=1.5s, **d6**=65 ms

PROCESSING


9. Process the recorded data with **xfb**
10. The resulting 2D spectrum can baseline corrected by **abs1, abs2**



4. **getprsol**: Get probe and solvent dependent parameters (corresponding pulses and power levels for the current experiment)
5. HMBC experiments must be run without sample spinning
6. Tune and match the probehead (**atma**)

ACQUISITION

7. Set the appropriate **ns** and **ds** for the experiment
8. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).



^1H - ^1H Noesy

PRELIMINARY SET-UP


1. Run a **conventional ^1H spectrum**. Note the **o1** and **SW** values to optimize spectral widths in the corresponding 2D experiment
2. Create a new dataset (**edc**) (or read and copy a pre-existing experiment) and read the standard BRUKER parameter set (**rpar**) to record a 2D ^1H - ^1H NOESY spectrum with **rpar NOESYPSW all** (the pulse program **noesygp** can be visualized in the **PulsProg** section or with the **edcpul** command). **getprsol**: Get probe and solvent dependent parameters (corresponding pulses and power levels for the current experiment)
3. NOESY experiments must be run without sample spinning
4. Tune and match the probehead (**atma**)

ACQUISITION

5. Set the appropriate **ns** and **ds** for the experiment
6. Start acquisition by **rga** and then **zg** (the expected experimental time is displayed with the **expt** command).

PROCESSING

7. Process the recorded data with **xfb**
8. The resulting 2D spectrum can baseline corrected by **abs1**, **abs2**




SPECIFIC PARAMETERS (acqparms)

By default, the following parameters are set to:

2 td in F2 (1K-2K), **1 td** in F1 (128w-256w)

ns= 2, **ds=16**

d1 =2s, **d8** (NOESY mixing time)= 0.1-0.2s (for large molecules and 0.4-0.5s for small ones)



Bibliography

<https://nmr.ucdavis.edu/sites/g/files/dgvnsk4156/files/docs/ucdavis-topspin3.2-userguide-october2015.pdf>

1D and 2D Experiments Step-by-Step Tutorial. Basic Experiments User Guide (Version 004) Bruker

1D and 2D Experiments Step-by-Step Tutorial Advanced Experiments User Guide (Version 002) Bruker

<http://www2.chem.uic.edu/nmr/downloads/BASHCD10/pdf/b4472.pdf>