

Time (s)	Calorimeter (oC)	Calorimeter - after cooling (oC)	Cu (oC)	Fe (oC)	Al (oC)
0					
30					
60					
90					
120					
150					
180					
210					
240					
270					
300					
330					
360					
365					
370					
375					
380					
385					
390					
395					
400					
405					
410					
415					
420					
435					
450					
465					
480					
495					
510					
525					
540					
555					
570					
585					
600					
630					
660					
690					
720					
750					
780					
810					
840					
870					

900					
930					
960					
990					
1020					
1050					
1080					